

KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody
Rabbit monoclonal antibody
Catalog # AGI2410**Specification****KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody - Product Information**

Application	WB, FC
Primary Accession	P17301
Reactivity	Human
Clonality	Monoclonal
Isotype	Rabbit IgG
Calculated MW	Predicted, 129 kDa, observed, 150 kDa
Gene Name	ITGA2
Aliases	ITGA2; Integrin Subunit Alpha 2; GPIa; VLAA2; HPA-5; CD49B; Integrin, Alpha 2 (CD49B, Alpha 2 Subunit Of VLA-2 Receptor); Very Late Activation Receptor Subunit Alpha-2; Alpha 2 Subunit Of VLA-2 Receptor; CD49 Antigen-Like Family Member B; Platelet Membrane Glycoprotein Ia; Human Platelet Alloantigen 5 2 3; Alpha 2 Integrin; Integrin Alpha-2; Glycoprotein Ia; Very Late Activation Protein 2 Receptor, Alpha-2 Subunit; Very Late Activation Receptor Alpha-2 Subunit; Human Platelet Alloantigen System 5; Platelet Glycoprotein GPIa; Platelet Antigen Br; VLA-2 Subunit Alpha; CD49b Antigen; VLA-2; BR
Immunogen	A synthesized peptide derived from human Integrin alpha 2

KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody - Additional Information

Gene ID	3673
Other Names	
Integrin alpha-2, CD49 antigen-like family member B, Collagen receptor, Platelet membrane glycoprotein Ia, GPIa, VLA-2 subunit alpha, CD49b, ITGA2, CD49B	

KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody - Protein Information**Name** ITGA2

Synonyms CD49B

Function

Integrin alpha-2/beta-1 is a receptor for laminin, collagen, collagen C-propeptides, fibronectin and E-cadherin. It recognizes the proline-hydroxylated sequence G-F-P-G-E-R in collagen. It is responsible for adhesion of platelets and other cells to collagens, modulation of collagen and collagenase gene expression, force generation and organization of newly synthesized extracellular matrix. (Microbial infection) Integrin ITGA2:ITGB1 acts as a receptor for Human echoviruses 1 and 8.

Cellular Location


Membrane; Single-pass type I membrane protein.

KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody - Protocols


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)


KO Validated Anti-Integrin Subunit Alpha 2 Rabbit Monoclonal Antibody - Images



Western blotting analysis using anti-Integrin subunit alpha 2 antibody (Cat#AGI2410). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Integrin subunit alpha 2 antibody (Cat#AGI2410, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-integrin subunit alpha 2 antibody (Cat#AGI2410). Integrin subunit alpha 2 expression in wild type (WT) and integrin subunit alpha 2 (ITGA2) knockout (KO) HSHC cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-integrin subunit alpha 2 antibody (Cat#AGI2410, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of Integrin subunit alpha 2 expression in HepG2 cells using anti-Integrin subunit alpha 2 antibody (Cat# AGI2410, 1:2,000). Green, isotype control; red, Integrin subunit alpha 2.