

Phospho-Ammonium transporter(T465) Antibody
Affinity Purified Rabbit Polyclonal Antibody (Pab)
Catalog # AP3770a**Specification**

Phospho-Ammonium transporter(T465) Antibody - Product Information

Application	DB,E
Primary Accession	E2CWJ2
Host	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Calculated MW	53698

Phospho-Ammonium transporter(T465) Antibody - Additional Information**Other Names**

AMT1;1 ammonium/proton co-transporter

Target/Specificity

This Phospho-Ammonium transporter-T465 Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide surrounding T465 of French bean Ammonium transporter.

Dilution

DB~~1:500

E~~Use at an assay dependent concentration.

Format

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Precautions

Phospho-Ammonium transporter(T465) Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Phospho-Ammonium transporter(T465) Antibody - Protein Information

Name E2CWJ2

Cellular Location


Cell membrane {ECO:0000256|RuleBase:RU362002}; Multi-pass membrane protein {ECO:0000256|RuleBase:RU362002}. Membrane {ECO:0000256|ARBA:ARBA00004141}; Multi-pass membrane protein {ECO:0000256|ARBA:ARBA00004141}

Phospho-Ammonium transporter(T465) Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Phospho-Ammonium transporter(T465) Antibody - Images



Dot blot analysis of Phospho-Ammonium transporter-T465 Antibody Phospho-specific Pab (Cat. #AP3770a) on nitrocellulose membrane. 50ng of Phospho-peptide or Non Phospho-peptide per dot were adsorbed. Antibody working concentrations are 0.6ug per ml.