

**Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL)
Recombinant Antibody
Catalog # APR10997****Specification**

Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) - Product Information


Application	FC, Kinetics, Animal Model
Primary Accession	NA
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1
Calculated MW	150 KDa

Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) - Additional Information**Target/Specificity**
oxLDL**Endotoxin**
< 0.001EU/ µg,determined by LAL method.**Conjugation**
Unconjugated**Expression system**
CHO Cell**Format**
Purified monoclonal antibody supplied in PBS, pH6.0, without preservative.This antibody is purified through a protein A column.**Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) - Protein Information****Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) - Protocols**


Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) - Images



Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%



The purity of Anti-oxLDL Reference Antibody (Okayama U. patent anti-oxLDL) is more than 95%, determined by SEC-HPLC.