

## Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody

Goat Polyclonal, Peroxidase (Horseradish) Catalog # ASR1884

### Specification

# Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Product Information

Description Host Conjugate Target Species Reactivity Clonality Application Application Note	Anti-CHICKEN IgG (H&L) (GOAT) Antibody Peroxidase Conjugated Goat Peroxidase (Horseradish) Chicken Chicken Polyclonal WB, E, IC ELISA 1:10,000-1:150,000;Western Blot 1:1,000-1:5,000;Immunochemistry 1:1,000-1:10,000
Physical State Host Isotype Target Isotype Buffer Immunogen	Lyophilized IgG IgG (H&L) 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 Anti-Chicken IgG (H&L) was produced by repeated immunization with Chicken IgG whole molecule in goat.
Reconstitution Volume Reconstitution Buffer Stabilizer Preservative	1.0 mL Restore with deionized water (or equivalent) 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free 0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!

## Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Additional Information

Shipping Condition Ambient

#### Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Chicken IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum, Chicken IgG and Chicken Serum.

### Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.



Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protein Information

## Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- <u>Western Blot</u>
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- <u>Cell Culture</u>

## Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Images

## Anti-CHICKEN IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Background

Anti-Chicken IgG antibody generated in goat detects specifically chicken IgG. This secondary antibody anti-Chicken is ideal for investigators who routinely perform ELISA, Sandwich ELISA, titration assays, western-blot, immunoprecipitation and more generally immunoassays.