

**Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody  
Catalog # ASR1961****Specification****Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody -  
Product Information**

Description	<b>Anti-GOLDEN SYRIAN &amp; ARMENIAN HAMSTER IgG (H&amp;L) (GOAT) (Min X MOUSE and RAT Serum Proteins)</b>
Host	<b>Goat</b>
Conjugate	<b>Unconjugated</b>
Target Species	<b>Armenian and Golden Syrian Hamster</b>
Application	<b>WB, IHC, E</b>
Application Note	<b>ELISA 1:20,000-1:100,000;Western Blot 1:2,000-1:10,000;Immunohistochemistry 1:1,000-1:5,000</b>
Physical State	<b>Liquid (sterile filtered)</b>
Host Isotype	<b>IgG</b>
Target Isotype	<b>IgG (H&amp;L)</b>
Buffer	<b>0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2</b>
Immunogen	<b>Armenian and Golden Syrian Hamster IgG, whole molecule</b>
Stabilizer	<b>None</b>
Preservative	<b>0.01% (w/v) Sodium Azide</b>

**Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody -  
Additional Information****Shipping Condition**

Wet Ice

**Purity**

Secondary Antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Golden Syrian and Armenian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Golden Syrian Hamster IgG and Armenian Hamster IgG. No reaction was observed against Mouse or Rat serum proteins. Specificity was confirmed using ELISA at less than 1% cross reactivity.

**Storage Condition**

Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.

**Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

**Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody -**

## **Protein Information**

### **Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)

### **Anti-Golden Syrian & Armenian Hamster IgG (H&L) Pre-Adsorbed Secondary Antibody - Images**