

### Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody

Goat Polyclonal, Peroxidase (Horseradish)
Catalog # ASR2455

#### **Specification**

## Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Product Information

Description Anti-GOLDEN SYRIAN HAMSTER IgG (H&L)

(GOAT) Antibody Peroxidase Conjugated

Host Goat

Conjugate Peroxidase (Horseradish)
Target Species Golden Syrian Hamster
Reactivity Hamster

Clonality Polyclonal Application WB, E, IC

Application Note ELISA 1:7,000-1:35,000;Western Blot

1:1,000-1:5,000;Immunochemistry 1:500-1:2,500

Physical State Lyophilized

Host Isotype IgG

Target Isotype IgG (H&L)

Buffer 0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Immunogen

Anti-Golden Syrian Hamster IgG whole molecule was produced by repeated

immunization with Golden Syrian Hamster

IgG whole molecule in goat.

Reconstitution Volume 1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Gentamicin Sulfate. Do NOT

add Sodium Azide!

# Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Additional Information

#### **Shipping Condition**

**Ambient** 

#### **Purity**

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Goat Serum, Golden Syrian Hamster IgG and Golden Syrian Hamster Serum. Diminished reactivity will occur against Armenian Hamster IgG.

#### **Storage Condition**

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C



or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

#### **Precautions Note**

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protein Information

Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Images

Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Background

Anti-Golden Syrian Hamster IgG whole molecule antibody generated in goat detects specificallyGolden Syrian Hamster IgG whole molecule. This secondary antibody anti-Golden Syrian Hamster is ideal for investigators who routinely perform ELISA, Sandwich ELISA, titration assays, western-blot, immunoprecipitation and more generally immunoassays.