

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody Goat Polyclonal, Peroxidase (Horseradish)
Catalog # ASR2882

Specification

Host

Conjugate

Clonality

Application

Target Species

Application Note

Physical State

Target Isotype

Host Isotype

Immunogen

Stabilizer

Preservative

Reconstitution Volume

Reconstitution Buffer

Buffer

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Product Information

Description F(ab) Anti-GS HAMSTER IgG (H&L) (GOAT)

Antibody Peroxidase Conjugated

Goat

Peroxidase (Horseradish) Golden Syrian Hamster

Polyclonal WB, E, IC

ELISA 1:10,000-1:50,000; Western Blot 1:1,000-1:10,000; Immunochemistry

1:500-1:2,500 Lyophilized IgG F(ab) IgG (H&L)

0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Golden Syrian Hamster IgG whole

molecule 1.0 mL

Restore with deionized water (or

equivalent)

10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

0.01% (w/v) Gentamicin Sulfate. Do NOT

add Sodium Azide!

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Additional Information

Shipping Condition Ambient

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Golden Syrian Hamster IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, papain digestion and chromatographic separation. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase and anti-Goat Serum. No reaction was observed against anti-Papain or anti-Goat IgG F(c).

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.



Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protein Information

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Fab Anti-Golden Syrian Hamster IgG (H&L) (Peroxidase Conjugated) Secondary Antibody - Images