

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody

Goat Polyclonal, Peroxidase (Horseradish)
Catalog # ASR2925

Specification

Physical State

Host Isotype

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Product Information

Description Anti-MOUSE IgM (H&L) (GOAT) Antibody

Peroxidase Conjugated

Host Goat

Conjugate Peroxidase (Horseradish)

Target Species
Clonality
Application

Mouse
Polyclonal
WB, E, IC

Application Note ELISA 1:10,000-1:50,000; Western Blot

1:1,000-1:10,000;Immunochemistry

1:500-1:2,500 Lyophilized

lgG

Target Isotype IgM (H&L)

Buffer 0.01 M Sodium Phosphate, 0.15 M Sodium

Chloride, pH 7.2

Immunogen Anti-Mouse IgM secondary antibody was

produced by repeated immunization with Mouse IgM whole molecule in rabbit.

Reconstitution Volume 2.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Gentamicin Sulfate. Do NOT

add Sodium Azide!

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Additional Information

Shipping Condition Ambient

Purity

Anti-Mouse IgM secondary antibody is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase and anti-Goat Serum. This reagent primarily recognizes Mouse IgM however some light chain cross reactivity will occur against other mouse immunoglobulins.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after



standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Protein Information

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Images

Anti-Mouse IgM (H&L) (Peroxidase Conjugated) Secondary Antibody - Background

Anti-Mouse IgM antibody generated in rabbit detects specifically Mouse IgM with some light chain cross reactivity against Mouse IgG. This secondary antibody anti-Mouse is ideal for investigators who routinely perform ELISA, Sandwich ELISA, titration assays, western-blot, immunoprecipitation and more generally immunoassays.